

Acquisition Core Project 3: Radiofrequency (RF) Pulse Development and Design for High Field MRI

Hui Liu^{1,2} and Gerald B. Matson^{1,3}

¹ Center for Imaging of Neurodegenerative Diseases, San Francisco, CA,
Northern California institute for Research and Education, and University of
California, San Francisco, CA

AIMS:

1. Develop RF pulses with immunity to B_1 inhomogeneity and to resonance offset
Initial focus on broadband (BB), shallow tip pulses for 3D MP RAGE
2. Expand simulations of magnetization transfer (MT) effects created by the new pulses
Initial focus on 3D MP RAGE simulations




BB Shallow Tip Pulses and 3D MP RAGE Simulations

- Pulse Generation
 - Development of optimal control routine for MatPulse patterned after that described by Skinner and Glaser in articles in JMR
 - Initial pulse (Parent pulse) submitted to optimal control routine in MatPulse
- Magnetization Transfer Simulations
 - Simulations based on formulation by Helms and Hagberg in Concepts in Magnetic Resonance




Parent Pulse Cascade – Null Pulse Version

First and Second Pulses



First pulse of parent pulse




First and second pulses of parent pulse


Trace diagrams of magnetizations experiencing different B_1 fields (color) and different resonance offsets (gray).

Parent Pulse Cascade – Null Pulse Version

Third and Fourth Pulses



First three pulses of parent pulse



All four pulses of parent pulse

The trace diagram shows that magnetizations experiencing different B_1 fields (color) are coalesced half way through the pulse, and again at the end of the pulse, while magnetizations experiencing different offsets (gray) are coalesced at the end of the pulse.

Optimal Control Menu in MatPulse

Op_Con_Non

BROADBAND B1-IMMUNE PULSE DESIGN by OPTIMAL CONTROL

Load B1 parameters Dwell (us): 3.7 Length (ms): 1.48 Pulse Points: 400

Optimization Parameters:

Pulse Type: Excite Tip (degrees): 10.0 ☒ Initiate ☐ Continue CALC

☐ Coalesced ☒ Non-Coalesced ☐ Linear: Factor Excitation (kHz): 0.2

Excitation Band Points: 41 B1 Immunity Range (+/-): 20 % Steps: 5

Limit B1 Maximum To: 7.0 uT ☐ Limit SAR SAR Factor (0.1 - 1.0): 0.5

Initial Step Size (0.1 - 1.0): 0.5 ☒ Average ☐ Fix Size ☐ Compute

Error Increase tolerance (1e-6 - 1e-3): 1e-4

Limit optimization by:

☐ Max iteration number: 64

☒ Residual error (%): 0.2

☐ Differential error (%%): 1e-4

☐ Halt if error increasing


STOP CLEAR

PB


The MatPulse menu showing the inputs for pulse optimization by optimal control



Example of Broadband (Non-Slice Selective) 10 Degree Tip RF Pulse Suitable for MP RAGE




Real and imaginary components of RF pulse.



Contour plot of pulse performance.

The pulse performs well over a B_1 range of $\pm 20\%$, and over a resonance offset range of ± 100 Hz. This range is suitable for 4.0 T.


Trace Diagram of the 10 Degree Pulse




Trace diagram of the 10 degree pulse (viewing the unit circle from above the z-axis) .

The trace diagram shows that magnetizations experiencing different B_1 field strengths (color) are coalesced halfway through the pulse, and again at the end of the pulse, while magnetizations at different resonance offsets (gray) are coalesced at the end of the pulse.

Example of Broadband (Non-Slice Selective) 30 Degree Tip RF Pulse Suitable for MP RAGE



Real and imaginary components of RF pulse.



Contour plot of pulse performance.

The pulse performs well over a B_1 range of $\pm 30\%$, and over a resonance offset range of ± 100 Hz. This range is suitable for 7.0 T.

Magnetization Transfer Menu for 3D MP RAGE in MatMRI

mr_parameters

PARAMETERS

EXECUTE

Tissue Parameters				Mag Trans Parameters			
Tissu	% H2O	T1 (ms)	T2 (ms)	T2* (ms)	% Bounc	T2 (us)	Ex Rate
Gray	80	1724	70	20	8	17	38.4
White	71	1043	65	20	16	16	43.3
CSF	100	4500	120	40	0	0	0
Other	0	0	0	0	0	0	0

Note: Unused parameters maybe left at zero

IMPORT Parameters from file

SAVE Parameters to file

RESET

CLEAR

mr_3dmprage

3D MP RAGE

Constant Tip Angle

EXECUTE

Inputs Parameters

Import from File

Inversion

Pulse Duration (us) / Strength (Hz) 1000 500

☒ Adiabatic (perfect)

☐ Imperfect: Tip Angle (deg)

☐ Use B1 Inhomogeneity (%)

Tissues (use Ctrl for multiple selections)

Gray

White

CSF

Other

RELOT

Excitation

Pulse Duration (us) / Saturation (delb) 1000 0.98

Tip Angle (deg) 10

Number of Pulses 22

Interpulse Delay (ms) 40

k-Space Center Pulse 0

Plot (multiple choices allowed)

☒ Plot Mz

☒ Plot Mxy

Options

☒ Include MT

Timings

TI (ms) 1200 Acquisition Period (ms) 880

TR (ms) 2400 Postacq Recovery (ms) 320


Save Export Parameters to File

CLEAR


The MatMRI menus showing the inputs for the magnetization transfer simulations



Simulations of an MP RAGE Experiment with and without Magnetization Transfer (MT) Effects



Simulation of an MP RAGE experiment showing Mz without MT effects



Simulation of an MP RAGE experiment showing Mz with MT effects

- MP RAGE experiments will have to be re-designed to accommodate the greater SAR and larger MT effects associated with the new shallow tip pulses.

Summary

- Broadband (non-slice selective) pulses are demonstrated with immunity to B_1 inhomogeneity and to resonance offset
- The pulses can be used for 3D MRI experiments such as 3D MPRAGE
- The new pulses are much longer than the rectangular pulses they replace, thus producing more SAR and inadvertent MT effects
 - MRI sequences will have to be re-designed with fewer pulses to accommodate the increased SAR and to mitigate the large MT effects
 - More k-space data must be collected following each pulse to avoid lengthening the duration of the MRI experiment
 - Corrections to the k-space data will have to be applied due to the long readouts

